



Genotyping of NDRG1 KO mice

1. Method used in our animal resource bank

This section describes our method we are routinely using for genotyping NDRG1 KO mice. This method is a little bit modified from the original (see the next section).

1.1) Primer sequences:

- Primer1: NDRG5'
 - Sequence: 5' - CCGCCTCTGTCAAATTAGTAGCTG -3' (24-mer)
- Primer2: NDRG3'
 - Sequence: 5' - GGGAGAGCTGAAGGCTGTTCTAGG -3' (24-mer)
- Primer3: NULL5'
 - Sequence: 5' - AGCAGGCTCTTAAAGCGGCTCC -3' (22-mer)

1.2) Reaction mixture:

| | For wild allele | For KO allele |
|------------------------------------|-----------------|---------------|
| | Tube 1 (µL) | Tube 2 (µL) |
| Water | 8 | 8 |
| Primer1 (NDRG5', 10 µM) | 0.5 | - |
| Primer2 (NDRG3', 10 µM) | 0.5 | 0.5 |
| Primer3 (NULL5', 10 µM) | - | 0.5 |
| Taq polymerase (U/µL) | 10 | 10 |
| DNA extracted from tail (purified) | 1 | 1 |
| total | 20 | 20 |

Taq polymerase: HotStarTaq Master Mix Kit (Qiagen). The enzyme is a chemically modified Taq polymerase for hot start PCR and needs 15-min incubation at 95 °C for activation. Master Mix contains enzyme, dNTP, Mg, etc at 2 x concentration. Please see Qiagen's website for details (<http://www1.qiagen.com/Products/Pcr/HotStarTaqSystem/HotStarTaqMasterMix.aspx>).

1.3) Thermal cycles:

| | | |
|-------|--------|--------------------------------------|
| 95 °C | 15 min | Enzyme activation and first denature |
| 94 °C | 15 sec | 35 cycles |
| 60 °C | 15 sec | |
| 72 °C | 20 sec | |
| 72 °C | 3 min | |
| 4 °C | ∞ | |

Thermal cycler: Veriti with 0.2mL tubes.

1.4) Product size:

Primers 1 and 2: approx. 200 bp for wild-type alleles

Primers 2 and 3: approx. 270 bp for KO alleles




Fig.1. Electropherogram of PCR products from heterozygous mice.

1.5) Reference

Okuda T, Higashi Y, Kokame K, Tanaka C, Kondoh H, Miyata T. Ndr1-deficient mice exhibit a progressive demyelinating disorder of peripheral nerves. *Mol Cell Biol.* 2004 May;24(9):3949-56. PMID: 15082788

Nbio_ID=062;05/30/2012